

Carpet pH Rinse

HILLYARD
ARSENAL®
35

Removes alkaline build up from carpet fibers and neutralizes them. Leaves a hard crystal residue that can be easily vacuumed out of the carpet. Low pH ensures removal of tough alkaline soils and some acid spills, such as coffee. Specially developed to be a part of a two-step system. 1: Pre-spray carpet with Hillyard Carpet Pre-Spray, HD Extraction or Clean Action II according to label directions. 2: Rinse with a solution of Carpet pH Rinse by applying through extraction machine.



Features & Benefits

- Removes and neutralizes alkaline buildup from carpet fibers.
- Neutralizes pH to prevent browning.
- Brightens the color of the carpet.

Directions

Arsenal II

1. For use with Hillyard Arsenal II Push Button Dispensers. Remove shipping cap and attach adapter cap to the one gallon of concentrate. Place into either a wire rack or cabinet.
2. Dispense diluted product into an extractor machine.

Arsenal Junior & Mobile

1. For use with Hillyard Arsenal Push Button Dispensers. Remove shipping cap and insert container into bottom of dispenser until it locks into place.
2. Dispense diluted product into a carpet extractor.

Arsenal Hil-Pac

Empty contents of pouch into 2 gallons of water.

Arsenal One

1. For use with the Hillyard Arsenal One dispensing systems. Remove shipping cap and insert this container into the dispensing system.
2. Secure dispensing system cap to refill container. Dispense product by selecting the product on the dispenser control panel and press the button to dispense product into bottle.
3. For bucket fill, attach bucket fill hose BEFORE dispensing product.

NOTICE: Over-wetting of carpet may cause shrinkage and browning.

Safety

See material safety data sheet and product label for safety information, handling and proper use.

HMIS

| | Concentrate/RTU |
|--------------|-----------------|
| Health | 2/0 |
| Flammability | 0/0 |
| Reactivity | 0/0 |

Technical Specifications

| | |
|---------------------|-------------|
| Color | Colorless |
| Scent | Mild |
| Appearance | Clear |
| pH (concentrate) | 2.00 - 3.00 |
| Non-Volatile Matter | 4.00% |
| Dilution Rate | 1:64 |

Availability

| Item | Pack |
|------------|-----------------------------------|
| HIL0083506 | 4 - "F" Style 1 Gallon Containers |
| HIL0083522 | 6 - 1/2 Gallon Containers |
| HIL0083594 | 36 - 4 oz. Packets |
| HIL0083525 | 4 - 2.5 Liter Containers |

HILLYARD

P.O. Box 909
St. Joseph, MO 64502
1-800-365-1555
www.hillyard.com